

BOTTLES


2696 BOTTLES SOLUTION, PLAIN TOOLED NECK

Cat. Ref	Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm	Qty per case
2696/21	250	65	147	18	24
2696/24	500	80	172	23	12
2696/28	1000	105	200	28	8
2696/29	2000	130	265	33	4
2696/30	3000	150	281	44	4
2696/32	5000	190	320	44	2
2696/34	9500	187	476	54	1
2696/35	10000	220	420	54	1
2696/36	13250	238	445	54	1
2696/40	20000	300	535	54	1



2741 BOTTLES WEIGHING, WITH INTERCHANGEABLE STOPPER

ITEM Code	Capacity MI.	Approx OD mm.	Approx. Ht. mm	Qty per case
2741/05	5	20	40	24
2741/08	15	25	60	12
2741/09*	*20	50	35	12
2741/10	25	30	65	12
2741/12*	*40	60	40	12
2741/14	60	40	90	6



* Squat Form

2771 BOTTLES WASH, COMPLETE WITH INTERCHANGEABLE STOPPER

Cat. Ref	Capacity MI.	Stopper	Qty per case
2771/21	250	24/29	8
2771/24	500	24/29	8
2771/28	1000	29/32	8

