




The Transparent Specialist

4133

**CYLINDERS, GRADUATED, SINGLE METRIC SCALE, WITH POUR OUT ,
WITH HEXAGONAL BASE CLASS B**

CAT. REF.	CAP. ML.	GRADU- -ATION INTERVAL ML	TOLER- -ANCE ±ML	APPROX. HEIGHT MM	PACKING I/C/M	MRP / PC	Group
4133/05	5	0.1	0.1	110	1/8/-	200	C
4133/06	10	0.2	0.2	140	1/8/-	225	A
4133/10	25	0.5	0.5	170	1/12/-	270	A
4133/13	50	1.0	1.0	200	1/12/-	295	A
4133/15	100	1.0	1.0	260	1/12/72	335	A
4133/21	250	2.0	2.0	335	1/8/32	625	A
4133/24	500	5.0	5.0	390	1/8/32	825	A
4133/28	1000	10.0	10.0	470	1/0/25	1250	C
4133/29	2000	20.0	20.0	570	1/0/9	2400	C



DESICCATORS

4192

DESICCATORS WITH COVER, KNOB TOP

CAT. REF.	SIZE MM	GROUND FLANGE APPROX. I.D. MM	PACKING I/C/M	MRP / PC	Group
4192/01	100mm	105	0/1/-	4500	D
4192/02	150mm	154	1/2/12	4950	D
4192/03	200mm	202	0/1/6	5500	A
4192/04	250mm	260	-/1/-	8250	A
4192/05	300mm	300	-/1/-	13950	D
4192/06	500mm	500	-/1/-	45000	E




NEW

4220*

DESICCATORS VACUUM, WITH TABULATED COVER, STOPPER, WITH PTFE SPINDLE

CAT. REF.	SIZE. MM	GROUND FLANGE APPROX. ID MM	INTER-CHANGABLE JOINT SIZE	PACKING I/C/M	MRP / PC	Group
4220/01	100mm	105	29/32	0/1/-	9225	D
4220/02	150mm	154	29/32	1/2/12	13300	D
4220/03	200mm	202	29/32	0/1/6	14025	D
4220/04	250mm	260	29/32	-1/-	22500	D
4220/05	300mm	300	29/32	-1/-	25625	D
4220/06	500mm	500	34/35	-1/-	75000	E



*Vacuum type Desiccators, when ground surface is cleaned & greased, will hold a vacuum of 500 mm of Mercury (Hg) over a 24 hours period.

DISHES

4251

DISHES, CRYSTALLIZING

CAT. REF.	SIZE MM	APPROX. O.D. MM	APPROX. HEIGHT MM	QTY PER CASE	MRP/ PC	Group
4251/1708	100	100	50	16	330	D
4251/2109	150	150	75	12	575	D
4251/2310	190	190	100	2	970	D



4271

DISHES, CULTURE, PETRI

CAT. REF.	SIZE MM	APPROX. O.D. MM	APPROX. HEIGHT MM	PACKING I/C/M	MRP/ PC	Group
4271/0803	50	50	17	0/48/-	130	E
4271/1303	80	80	17	0/48/-	165	A
4271/1703	100	100	17	0/48/-	200	A
4271/2104	150	150	20	0/20/-	495	E
4271/2404	200	200	20	0/8/-	1075	E

