

8193

PIPETTES, SEROLOGICAL (WHITE MARKING)

CAT. REF.	CAP. ML.	GRADUATION INTERVAL ML	TOLERANCE ± ML	PACKING I/C/M	MRP/ PC	Group
8193/61*	0.1	0.01	0.01	1/8/-	110	E
8193/62*	0.2	0.01	0.01	1/8/-	105	E
8193/01	1.0	0.01	0.01	1/8/-	95	E
8193/03	2.0	0.02	0.02	1/8/-	110	E
8193/05	5.0	0.05	0.05	1/8/-	120	E
8193/06	10.0	0.1	0.1	1/8/-	125	E
8193/10	25.0	0.2	0.2	1/8/-	190	E

These pipettes are graduated to tip.



8211

PIPETTES, TRANSFER, VOLUMETRIC, CLASS A

CAT. REF.	CAP. ML.	TOLERANCE ± ML	PACKING I/C/M	MRP/ PC	Group
8211/01*	1*	0.007	1/8/-	220	E
8211/03*	2*	0.01	1/8/-	225	E
8211/05	5	0.015	1/8/-	250	E
8211/06	10	0.02	1/8/-	305	E
8211/09	20	0.03	1/8/-	390	E
8211/10	25	0.03	1/8/-	420	E
8211/13	50	0.05	1/8/-	480	E
8211/15	100	0.08	1/8/-	715	E

* not covered by IS specifications.



8213

PIPETTES, TRANSFER, VOLUMETRIC, CLASS B

CAT. REF.	CAP. ML.	TOLERANCE ± ML	PACKING I/C/M	MRP/PC	Group
8213/01	1	0.015	1/8/-	100	E
8213/03	2	0.02	1/8/-	100	E
8213/05	5	0.03	1/8/-	105	E
8213/06	10	0.04	1/8/-	125	E
8213/09	20	0.06	1/8/-	150	E
8213/10	25	0.06	1/8/-	155	E
8213/13	50	0.1	1/8/-	180	E
8213/15	100	0.16	1/8/-	280	E

